

Rabbit Anti-Mouse IgG1 Monoclonal Antibody

Cat. No.: C52760

Summary

Description	Rabbit anti-Mouse IgG1 is a secondary antibody.
Clonality	Monoclonal
Host Species	Rabbit
Target Species	Mouse
Isotype	IgG

Property

Purity	>95%
Purification	Affinity purified
Sterility	0.2 µM filtered
Storage	Store at 4°C within a week. For longer storage, aliquot and store at -20°C.

Applications

Application	WB, IP, FC, ELISA
--------------------	-------------------

Target

Target	Mouse IgG1
Alternative Name	Immunoglobulin Gamma 1; IgG1

For lab research use only, not for diagnostic, therapeutic or any *in vivo* human use.