

Recombinant Anti-Emb Antibody (V3S-0522-YC5905)

Cat. No.: V3S-0522-YC5905

Summary

Description	This product is a monoclonal antibody derived from Rat (<i>Rattus norvegicus</i>), which can specifically recognize Embigin. The antibody is expressed with mammalian cell transient expression system, serum-free and purified by affinity chromatography. The purity and integrity are tested via SDS-PAGE and SEC-HPLC analysis. Given an antigen, additional QC measures are also desired such as affinity testing and binding validation. Specifically, the antibody is provided in multiple formats for in vivo and in vitro assays. The <i>In vivo</i> version features greater than 95% purity, ultra-low endotoxin levels (<1 EU/mg or 0.1 EU/mg), and is preservative, stabilizer, and carrier protein-free.
Clonality	Monoclonal
Host Species	Rat
Target Species	Mouse
Immunogen	Mouse Embigin overexpressing cells
Isotype	IgG2b
Secondary Antibody	C56110; C40224; C74571; C22227; C56421; C32514; C47763; C63764; C41555; C10423; C53263; C44755

Property

Expression Species	HEK293F or CHO
Conjugation	None
Purity	>95%, determined by SDS-PAGE and/or SEC-HPLC
Endotoxin	<1 EU/mg, determined by LAL method
Purification	Protein A affinity purified
Sterility	0.2 µM filtered
Formulation	PBS, pH 7.4
Preservation	No preservatives
Stabilizer	No stabilizers
Storage	Store at 4°C within a week. For longer storage, aliquot and store at -20°C.

For lab research use only, not for diagnostic, therapeutic or any *in vivo* human use.

Applications

Application FC; WB

Application Notes The antibody is recommended for detection of Emb by FC, WB assays.

Target

Target Emb

Alternative Name Gp70; AL022799

Gene ID [13723](#)

UniProt [P21995](#)

Introduction Embigin is a transmembrane glycoprotein that is a member of the immunoglobulin superfamily. The encoded protein may be involved in cell growth and development by mediating interactions between the cell and extracellular matrix.

Research Area Immunology

Related Disease Multiple sclerosis, psoriasis, rheumatism, chronic noninfectious uveitis, and glomerulonephritis

For lab research use only, not for diagnostic, therapeutic or any *in vivo* human use.

SUITE 203, 17 Ramsey Road, Shirley, NY 11967, USA

Tel: 1-631-381-2994 Fax: 1-631-207-8356

Email: info@creative-biolabs.com

Page 2 / 2

© Creative Biolabs All Rights Reserved