

Anti-ALCAM Neutralizing Antibody (V3S-0622-YC726)

Cat. No.: V3S-0622-YC726

Summary

Description	This product is a human monoclonal antibody that reacts with ALCAM. The antibody is expressed with mammalian cell transient expression system, serum-free and purified by affinity chromatography. The purity and integrity are tested via SDS-PAGE and SEC-HPLC analysis. Given an antigen, additional QC measures are also desired such as affinity testing and binding validation. Specifically, the antibody is provided in multiple formats for in vivo and in vitro assays. The <i>In vivo</i> version features greater than 95% purity, ultra-low endotoxin levels (<1 EU/mg or 0.1 EU/mg), and is preservative, stabilizer, and carrier protein-free.
Clonality	Monoclonal
Host Species	Human
Target Species	Human
Isotype	IgG1
Isotype Control	C67002; C43061
Secondary Antibody	C32400; C75370; C10513; C51635; C45530

Property

Expression Species	HEK293F or CHO
Conjugation	None
Purity	>95%, determined by SDS-PAGE and/or SEC-HPLC
Endotoxin	<1 EU/mg, determined by LAL method
Purification	Protein A affinity purified
Sterility	0.2 µM filtered
Formulation	PBS, pH 7.4
Preservation	No preservatives
Stabilizer	No stabilizers
Storage	Store at 4°C within a week. For longer storage, aliquot and store at -20°C.

Applications

For lab research use only, not for diagnostic, therapeutic or any *in vivo* human use.

Application	FuncS
Application Notes	The antibody is recommended for detection of ALCAM by FuncS assay.

Target

Target	ALCAM
Alternative Name	ALCAM; activated leukocyte cell adhesion molecule; MEMD; CD16
Gene ID	214
UniProt	Q13740
Introduction	This protein is activated leukocyte cell adhesion molecule (ALCAM), also known as CD166 (cluster of differentiation 166), which is a member of a subfamily of immunoglobulin receptors with five immunoglobulin-like domains (VVC2C2C2) in the extracellular domain. This protein binds to T-cell differentiation antigen CD6, and is implicated in the processes of cell adhesion and migration. Multiple alternatively spliced transcript variants encoding different isoforms have been found.
Research Area	Immunology

For lab research use only, not for diagnostic, therapeutic or any *in vivo* human use.