

Recombinant Anti-TNFSF15 Antibody (V3S-1022-YC1225)

Cat. No.: V3S-1022-YC1225

Summary

Description	This product is a recombinant mouse antibody provided by CreativeBiolabs. The antibody is specific for TNF Superfamily Member 15, the protein encoded by this gene is a cytokine that belongs to the tumor necrosis factor (TNF) ligand family. It can be used for TNFSF15 detection in Enzyme-Linked Immunosorbent Assay. It was expressed in mammalian cells (293F or CHO) with antibody encoding genes and purified by affinity chromatography. Each lot of this antibody is quality control tested by SDS-PAGE and SEC-HPLC analysis. For highly sensitive assays, we recommend the ultrapurified form of the product, which has a lower endotoxin limit than standard antibody, less than 1 EU/mg or even 0.1 EU/mg.
Clonality	Monoclonal
Host Species	Mouse
Target Species	Human
Isotype	IgG
Isotype Control	C35500
Secondary Antibody	C37557; C41360; C32672; C10001

Property

Expression Species	HEK293F or CHO Cell line
Conjugation	Unconjugated, also available for Biotin, HRP, FITC and PE-labeled form.
Purity	>95%, determined by SDS-PAGE and SEC-HPLC
Endotoxin	<1 EU/mg, determined by LAL method
Purification	Purified with Protein A or G affinity chromatography
Sterility	0.2 µM filtered
Formulation	PBS, pH 7.4
Preservation	No preservatives
Stabilizer	No stabilizers
Storage	Store at 4°C within one or two weeks. Store at -20°C for long term. Avoid repeated freeze/thaw cycles. Refer to the COA file for specifics.

For lab research use only, not for diagnostic, therapeutic or any *in vivo* human use.

Applications

Application ELISA

Target

Target TNFSF15

Alternative Name TNF Superfamily Member 15; Vascular Endothelial Cell Growth Inhibitor; TNF Ligand-Related Molecule 1; Tumor Necrosis Factor (Ligand) Superfamily; Member 15; Vascular Endothelial Growth Inhibitor-192A; TNF Superfamily Ligand TL1A; VEGI; TL1

Gene ID [9966](#)

UniProt [O95150](#)

Introduction This protein is abundantly expressed in endothelial cells, but is not expressed in either B or T cells. The protein encoded by this gene is a cytokine that belongs to the tumor necrosis factor (TNF) ligand family.

Research Area Cardiovascular; Immunology

For lab research use only, not for diagnostic, therapeutic or any *in vivo* human use.